



COTTON FACE MASK TECHNICAL CHART - PPE



CODE	MS01
MATERIAL	100% DOUBLE TWISTED COTTON 165 gr.mq nr. 56 number of thread in warp and weft. Silicone coating RE: N 1100 allows 500 washings at 90°/194F
MANUFACTURER	PASTELLI S.R.L. – CAMPO SAN MARTINO – ITALIA
PACKAGING	10 PIECES NOT STERILISED – ASSORTED COLORS

DESCRIPTION
<ul style="list-style-type: none">• The mask is not a sterilised product.• The materials are tested and do not present any kind of toxicity.• The mask may undergo high temperature washing at 90°/194F.• The mask MS02 is a PPE able to fulfill the characteristics required by the WHO.• The mask has no deadline.

<ul style="list-style-type: none">• USE The cotton mask must be disinfected before the use.
INSTRUCTION FOR USE <ul style="list-style-type: none">• Washing the mask with high temperature water at 60°/140F and a specific additive for disinfection, before using.• Before using the mask, wash carefully the hands with soap and water, or with an alcoholic solution.• Covering mouth and nose with the mask, make sure of the intactness of it and that it fits snugly at the face.• Avoid touching the mask while you are wearing it.• Replace the mask with a new one after the use, touch only the elastics while you are changing the mask without getting in contact with the frontal part of the mask. Have care to put the used one in a bag before the washing process with disinfection.• You can reuse the mask only after having followed the correct washing system. We suggest to put the washed mask in a single-use bag.• The mask must not be used when the fabric is not intact anymore.• When the mask is used by an underage child, the application, the putting off and the correct use must be managed by an adult. <p style="text-align: center;"><u>The mask is for a personal use only; do not share it with other people</u></p>