

Sensorguide

Sensormärke	Artikellnummer	Färg	Size
Dentalmind Digital X-Ray II	C01	Violet	#2
Dentalmind Digital X-Ray II	C01	Cerise	#2
Dentalmind Digital X-Ray II	C01	Yellow	#2
Dentalmind Digital X-Ray I	C03	Green	#1
Dentalmind Digital X-Ray I	C03	Blue	#1
SuniRay	C01	Violet	#2
SuniRay	C01	Cerise	#2
SuniRay	C01	Yellow	#2
SuniRay	C03	Green	#1
SuniRay	C03	Blue	#1
Instrumentarium Snapshot	C01	Violet	#2
Instrumentarium Snapshot	C01	Cerise	#2
Instrumentarium Snapshot	C01	Yellow	#2
Instrumentarium Snapshot	C03	Green	#1
Instrumentarium Snapshot	C03	Blue	#1
Soredex Toto	C01	Violet	#2
Soredex Toto	C01	Cerise	#2
Soredex Toto	C01	Yellow	#2
Soredex Toto	C03	Green	#1
Soredex Toto	C03	Blue	#1
Schick CDR, Schick Elite	C01	Violet	#2
Schick CDR, Schick Elite	C01	Cerise	#2
Schick CDR, Schick Elite	C01	Yellow	#2
Schick CDR, Schick Elite	C03	Green	#1
Schick CDR, Schick Elite	C03	Blue	#1
Gendex GXS-700	C01	Violet	#2
Gendex GXS-701	C01	Cerise	#2
Gendex GXS-701	C01	Yellow	#2
Gendex GXS-702	C03	Green	#1
Gendex GXS-703	C03	Blue	#1
Mdx3 (digitales)	C01	Violet	#2
Mdx3 (digitales)	C01	Cerise	#2
Mdx3 (digitales)	C01	Yellow	#2
Mdx3 (digitales)	C03	Green	#1
Mdx3 (digitales)	C03	Blue	#1
Plandent Prosensor	C01	Violet	#2
Plandent Prosensor	C01	Cerise	#2
Plandent Prosensor	C01	Yellow	#2
Plandent Prosensor	C03	Green	#1
Plandent Prosensor	C03	Blue	#1

Digitala sensorer:

Violet: Liggande & Stående bitewing | **Size #2**

Cerise: Front & apikal | **Size #2**

Green: Liggande & Stående bitewing | **Size #1**

Blue: Front & apikal | **Size #1**

Yellow: Endo | **Size #2**

Film- / bildplattor

Violet / F01: Liggande & Stående bitewing.

Cerise / F01: Front & apikal